

FINISHED PRODUCT SPECIFICATION

PRODUCT NAME	ZEAXANTHIN 20.0% POWDER	SPEC:FPS/B01/1118-01
		Effective Date: 04-10-2021

Botanical/ Scientific name	Tagetes erecta	Plant part used	Petals of Marigold
Standardization	Zeaxanthin	Production Method	Drying
Shelf Life	3 years	Extraction Solvent	Ethanol

SL No	Protocol	Test Parameter	Test Limit			
Physic	Physical					
01	Organoleptic	Description	Orange to red color* powder with characteristic odor			
02	TM/1118	Identification	To Comply by HPLC			
03	USP	Solubility	Partially soluble in alcohol Insoluble in water			
04	USP <921>	Moisture content	Not more than 5.0%w/w			
05	USP <616>	Untapped bulk density	Between 0.40 to 0.70g/ml			
06	USP <616>	Tapped bulk density	Between 0.50 to 0.80g/ml			
07	USP <786>	Particle size	100% passing through 40 mesh			
Chemi	cal					
08	TM /1118	Content of Zeaxanthin by HPLC	Not less than 20.0% w/w			
Heavy	Heavy Metal					
	USP <233>Procedure 1	Lead	Not more than 3.0 ppm			
09	USP <233>Procedure 1	Arsenic	Not more than 1.0 ppm			
0)	USP <233>Procedure 1	Cadmium	Not more than 0.5 ppm			
	USP <233>Procedure 1	Mercury	Not more than 0.1 ppm			
Microl	Microbiological					
	USP <2021>	Total plate count	Not more than 3000cfu/g			
	USP <2021>	Total yeast and mold count	Not more than 100cfu/g			
	USP <2022>	Escherichia coli	Absent / 10g			
10	USP <2022>	Salmonella spp.	Absent / 10g			
	USP <2021>	Coliforms	Absent / 10g			
	USP <2022>	Staphylococcus aureus	Absent / 10g			
	USP <62>	Pseudomonas aeruginosa	Absent / 10g			























FINISHED PRODUCT SPECIFICATION

ADDITIONAL INFORMATION		
Sanitizing treatment	Non –irradiated and not treated with ETO	
GMO Status	GMO free	
Standard Packing	25 Kg in single HDPE container.	
Certifications	ISO 22000:2005, ISO 9001:2015, GMP, HALAL certified	
Storage condition	Store in a cool and dry place away from sunlight and heat. The preferred storage is below 30 ⁰ C.	

^{*} The product is of herbal origin. There may be minor color difference because of the geographical and seasonal variations of the plant material.

This is a computer generated copy. Hence signature not required



Where Quality Matters.....



















